



Instructions

Welcome to Peace of Mind

Walabot HOME is the first automatic fall detection device that doesn't need to be worn.

Just mount it on the wall, and Walabot HOME's revolutionary sensing technology will detect if you fall and will automatically call for help.

What's In the Box



Walabot HOME device



Instruction Manual




Power adapter



Double-sided adhesive strips

IMPORTANT: Watch our video installation guide to learn how to setup your Walabot HOME correctly



Helpful information you will learn in the video:

- ✓ Choosing the right location
- ✓ Mounting Walabot HOME securely
- ✓ Connecting with your caregiver
- ✓ Testing your device

Watch it now at

walabot.com/go

(We recommend you have your caregiver present during installation)



Need Assistance with Installation?

Contact us at:

 Email: [**homesupport@walabot.com**](mailto:homesupport@walabot.com)

 Or Call Us at: **1 (844) 868-1799**

Your safety and peace of mind matter to us. Whatever your question, we're here.

Warnings

Disclaimers

1. The Walabot HOME device must be installed in the room that you intend to monitor.
2. The device must be connected to the Walabot HOME App in order to communicate with the emergency contact.
3. The size of a room the Walabot HOME device can support is 10ft x 10ft.
4. The Walabot HOME device may not function properly if it is not installed in the correct location and in accordance with the instructions.
5. In order to give you time to respond to a false alarm, the device permits a short time delay between a fall and contacting the emergency contact. This delay can be up to two minutes.
6. Presence of pets or additional people in the room during a fall might cancel/ prevent the fall detection
7. Several people in the room at the same time can trigger a false alarm.
8. Do not install the device inside or over the bathtub.
9. The device will not work for anyone who is lying in a bathtub full of water. In this case, the device may trigger a false alarm.
10. Your Walabot HOME device must be paired to the emergency contact phone number you have selected.
11. Not all falls can be detected and the device shouldn't be used as a replacement to a medical device.
12. The device supports a 140 degrees field of view. A fall will not be detected outside of the field of view.
13. The device will not always work if there is no direct line of sight to the person who falls.
14. You should test the device and recommend repeating the test in different locations around the room to make sure the device covers different areas of the room.
15. Repeat the test fall once a month to make sure the device still functions correctly
16. The room being monitored's floor must have the same floor height to ensure falls are detected accurately.

Warnings

1. Use and handle your Walabot HOME device carefully.
2. The Walabot HOME device can break if dropped.
3. Using a power adaptor other than the one provided can cause malfunction or potentially damage the Device.
4. Opening or tampering with the Device voids the warranty.
5. Keep the device away from water. Water damage is not covered by the warranty.
6. The device will not function during a power outage.
7. The device will not transmit an alert to an emergency contact if the WiFi in your home is not working.
8. You must follow the installation and use instructions or the device will not work.
9. Do not mount on wallpaper because the device may fall off the wall or damage the wallpaper.

10. Follow the instructions carefully before removing the device from the wall.

FCC Regulations

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this Device may not cause harmful interference, and (2) this Device must accept any interference received, including interference that may cause undesired operation. Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital Device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.
- To maintain compliance with FCC's RF Exposure guidelines, this equipment should be installed and operated with a minimum distance of 20cm between the radiator and your body.

Use only the supplied antenna.

Changes or modifications made to this equipment not expressly approved by Vayyar Imaging Ltd. may void the FCC authorization to operate the equipment.

FCC ID: 2AHIS-WH21BBUS02 | Contains FCC ID: 2AHIS-VHOME.

IC Regulations

IC Caution: RSS-Gen Issue 4 December 2014"/"CNR-Gen 4e Décembre 2014:

- English:

This device complies with Industry Canada licence-exempt RSS standard(s).

Operation is subject to the following two conditions: (1) This device may not cause interference, and (2) This device must accept any interference, including interference that may cause undesired operation of the device.

- French:

Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes:

- (1) l'appareil ne doit pas produire de brouillage, et;
- (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'encourager le fonctionnement.

IC: 21498-WH21BBUS02 | Contains IC: 21498-VHOME.